

CERAGEM V6

THERAPEUTIC THERMAL MASSAGER @HOME

FDA Cleared Medical Claims

Intensive Massage Along **the Spinal Line to Toes**



Temporary relief of minor muscle and joint pain stiffness

Temporary relief of minor joint pain associated with arthritis

Temporary increase in local circulation where applied

Relaxation of muscles

Brand New



○ White

● Black

● Camel Brown



NEW

Intensive Close-Up Massage Along the Spine (Acupressure Massage)



PATENTED

Customized Massage Tailored to User's Spinal Length and Curvature



NEW

Promote Blood Circulation by Air Cell Massagers



UPGRADE

22 Massage Modes Developed by Well-Life Medical Research Center



NEW | PATENTED

Adjustable 3 Massage Speed Settings Fast, Normal, and Slow



NEW | PATENTED

Auto-Start Upon Lying Down



NEW

Easy Sliding Feature for Easier Folding and Unfolding

CERAGEM V6

Class II Medical Device with 4 FDA Cleared Medical Claims


Intensive Massage Along the Spinal Line to Toes

- + **Concentrated massage** at acupressure points!
- + The Heated Massage Rollers move upward to **Stretch out the spine!**
- + **Relieve tension** and loosen stiff muscles!
- + **Intensive thermal therapy** up to 149°F!
- + **Promote blood circulation** of legs by Air Cell Massagers!

Brand New



POINT 1.



Intensive Thermal and Concentrated Massage Therapy Helps Relieve Muscle Pain

The specialized heated massage rollers provide a concentrated and seamless massage at your acupressure points by automatically adjusting its vertical location. This is done by the first rollers moving in a flexible motion, as they are designed separately from the second.

POINT 3.

Customized Massage for Each User's Unique Body Line

Heated massage rollers scan and analyze user's spinal length and curvature to provide a customized decompressive massage for all kinds of body shapes from children to adults.

POINT 2.



The V6 Not Only Provides an Intensive Massage to the Spinal Region, but the Air Cell Massagers also Promote Circulation to the Legs

4 Air Cells provide an intensive massage to the legs.
- Repetitive and powerful pneumatic inflation and deflation
- Specialized 4 Air Cell sections and 3 massage modes equipped.

POINT 4.

22 Massage Modes Developed by Well-Life Medical Research Center

17 pre-programmed auto modes including the newly added Acupressure Mode, Golf Mode and others. Plus 5 custom modes including the Semi-auto mode and the Manual mode to suit your body condition.

EXPERIENCE THE NEW CERAGEM V6

Your Everyday Answer to Recharge



Upgraded Features of CERAGEM V6

- Heated Massager Heats up to 140°F
- 3-Level Massage Speed Settings
- 9-Level Massage Intensity Settings
- 20 Different Ceragem Sound Settings
- Easy Folding and Unfolding by Easy Sliding Feature
- Low Noise/ Anti-Vibration While in Usage
- Auto-Start
- Voice Guidance to Ease the Usage



* In V6 User Manual : Heated massage rollers (Internal projector), Heated massager (Abdominal vibration projector)

Website | www.ceragemus.com Customer Service | +1 (213) 480-7070 E-mail | info@ceragemus.com

CERAGEM